ADVOCATE®

Redi-Code⁺

Quick Reference Guide

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Dear ADVOCATE REDI-CODE Plus System Owner:

Thank you for choosing ADVOCATE Redi-CodePlus Blood Glucose Monitoring System to take control of your diabetes. This system is helpful in tracking your blood glucose values. Your Quick Reference Guide assists you to set up the system and helps you to run the test. But if you'd like more information and additional features, please check the User Manual. Please read this Guide carefully before use, or you can read the User Manual for complete information.

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
Customer Service/Servicio al Cliente:

9:00 am - 5:00pm ET
Monday - Friday

At all other times, you should contact your health care professional for assistance.

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Test Strips
Tira Reactivas

Control Solution
Solución de Control

1 INSERT A TEST STRIP TO TURN ON THE METER.

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2 CHOOSE TEST MODES. When is flashing press M button to choose test modes (pre-meal/post-meal/general)

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3 While is flashing on the display, apply blood to the absorbent hole of the test strip until the confirmation window is completely covered.

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The blood should be completely filled in the confirmation window. If you find that the confirmation window is not filled with blood when the meter is counting, NEVER try to add more blood to the test strip. Discard the test strip and retest with a new one.

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4 After the meter counts to 5, your test result appears along with Date/Time/Result indicator and is automatically stored in the memory.

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Diabetic Supply of Suncoast, Inc.