# ADVOCATE<sup>®</sup>

*Redi-Code<sup>+</sup>*

***Owner's Manual***      ►►►  
***Manual del Usuario***

Rev. 10/2011

# Table of Contents

<b>Important Safety Instructions</b>	03
<b>Alternate Site Testing (AST)</b>	06
<b>Important Information about Your System</b>	08
Intended Use	08
Test Principle	08
<b>Introduction of Your System</b>	10
System Contents	10
Meter Components	11
Meter Display	12
Test Strip	13
<b>Meter Set-Up</b>	14
Install Batteries	14
Set the Meter	16
Use the Alarm	21
Cancel the Alarm	21
<b>Quality Control</b>	22
Important Control Solution Test Information	22
Run a Control Solution Test	23
<b>Blood Glucose Test</b>	26
Prepare for the Test	26
Run a Blood Glucose Test	28
Expected Blood Glucose Values	31
<b>Limitations</b>	32
<b>Comparison of Your Meter and Laboratory Results</b>	33

# Table of Contents

<b>Meter Memory</b> .....	35
7-, 14-, 28-, 60- Day Average Results .....	35
Test Results .....	37
<b>Your System Maintenance</b> .....	39
Cleaning and Disinfection .....	39
Storage .....	48
<b>Messages and Troubleshooting</b> .....	49
<b>Specifications</b> .....	54
<b>Statements</b> .....	55
<b>Warranty</b> .....	56

### Important Safety Instructions

- The meter and lancing device are for single patient use. Do not share them with anyone including other family members! Do not use on multiple patients! <sup>\*1,2</sup>
- All parts of the kit are considered biohazardous and can potentially transmit infectious diseases, even after you have performed cleaning and disinfection. <sup>\*1,2</sup> Details of the cleaning and disinfecting procedure could be found in the section Cleaning and Disinfection in this manual.
- Due to the potential biohazardous concern, you should wash hands thoroughly with soap and water after handling the meter, lancing device, or test strips.
- Severe dehydration and excessive water loss may cause false low results. If you believe you are suffering from severe dehydration, consult your healthcare professional immediately.
- If your blood glucose results lower or higher than usual, and you do not have symptoms of illness, first repeat the test. If you have symptoms or continue to get results higher or lower than usual, follow the treatment advice of your healthcare professional.
- If you are experiencing symptoms that are not consistent with your blood glucose test results and you have followed all instructions described in this Owner's Manual, call your healthcare professional.
- Inaccurate results may occur in severely hypotensive individuals or patients in shock. Readings which are lower than actual values may occur for individuals experiencing a hypoxia state, or a hyperglycemic hyperosmolar state, with or without ketosis. Critically ill patients should not be tested with blood glucose meters.

# Important Safety Instructions

- Do not use this meter in a dry environment, especially if synthetic materials are present. Synthetic clothes, carpets, etc., may cause damaging static discharges in a dry environment. Avoid touching the insertion slot for the test stick if in a dry environment.
- Do not use this meter near cellular or cordless telephones, walkie talkies, garage door openers, radio transmitters, or other electronical equipment that are sources of electromagnetic radiation, as these may interfere with the proper operation of the meter.

## ■ Lipemic Effects:

Cholesterol up to 500 mg/dL or triglycerides up to 2000 mg/dL do not significantly affect the results. However, the glucose test may be affected beyond these levels.

## ■ Interferences:

When tested following NCCLS guidelines<sup>\*3</sup>, bilirubin, creatinine, methyl dopa, galactose, maltose, xylose, ibuprofen, and salicylate at normal therapeutic concentrations do not significantly affect glucose results. However, these levels of the following interferents in blood may cause inaccurate test results:


- Acetaminophen  $\geq 12.5$  mg/dL (Therapeutic level is 1.2-3.6 mg/dL)
- Ascorbic acid  $\geq 7.5$  mg/dL (Therapeutic level is 0.4-2.1 mg/dL)
- Dopamine  $\geq 3$  mg/dL (Therapeutic level is 0.04 mg/dL)
- L-dopa  $\geq 4$  mg/dL (Therapeutic level is 0.02-0.3 mg/dL)
- Tolbutamide  $\geq 150$  mg/dL (Therapeutic level is 3.6-7.2 mg/dL)
- Uric acid  $\geq 20$  mg/dL (Normal level is 2.6-7.2 mg/dL)<sup>\*4</sup>
- Gentisic acid  $\geq 25$  mg/dL (Therapeutic level is 0.2-0.6 mg/dL)
- Tolazamide  $\geq 15$  mg/dL (Therapeutic level is 2.0-2.5 mg/dL)<sup>\*5</sup>
- Mannose  $\geq 250$  mg/dL (Therapeutic level is 1.15 mg/dL)
- Ibuprofen  $\geq 50$  mg/dL (Therapeutic level is 1.0-7.0 mg/dL)

## References

1. FDA Public Health Notification: Use of Fingertick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens: Initial Communication (2010) <http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/ucm224025.htm>
2. CDC Clinical Reminder: Use of Fingertick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens(2010) <http://www.cdc.gov/injectionsafety/Fingertick-DevicesBGM.html>
3. NCCLS, Vol. 25, No. 27, 2005, EP7-A2, Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition.
4. Tietz Textbook of Clinical Chemistry and Molecular Diagnosis (4th Ed.), Burtis, Ashwood & Bruns (Eds), Elsevier Saunders, 2005; 2301.
5. Bioavailability of Tolazamide from Tablets: Comparison of *In Vitro* and *In Vivo* Results, P. G. Welling, R. B. Patel, U. R. Patel, W. R. Gillespie, W. A. Craig, K. S. Albert, 2006.

## Alternate Site Testing (AST)

This system gives you to test on the palm, the forearm, the upper arm, the calf, or the thigh with the equivalent results to fingertip testing. These areas of your body are called alternate sites other than fingertip. Since nerve endings (receptors) are not so condensed at alternate sites, tests performed there will be less painful.

### ***When you should use AST***

Alternate site results may be different from fingertip results when glucose levels are changing rapidly. Therefore, **use AST only when it has been more than two hours after any one of the following:**

- A meal.
- Diabetes medication.
- Exercise.

When testing blood glucose during or immediately after meals, medication, physical exercise, or any other event, **take blood sample from your finger only.**

### ***Do NOT use blood samples from AST when***

- You think your blood glucose is low.
- You have been diagnosed with hypoglycemic unawareness.
- The AST results do not agree with the way you feel.
- You are ill or stressed.
- It is within two hours of eating a meal, taking insulin or exercising.

**AST results serve as references only. Consult your healthcare professional before performing alternate site testing (AST).**

## ***Before you get a drop of blood for AST***

- Rub the puncture site for 20 seconds before penetration for increasing the accuracy.

## ***Important Notices***

- Alternate site testing can be used only during steady-state blood glucose conditions.
- Measurements from alternative site testing should never be used to calibrate a continuous glucose monitor (CGM) or entered into insulin dose calculators for insulin dosing recommendations.

## ***Important Information about Your System***

### ***Intended Use***

Advocate® Redi-Code<sup>+</sup> Blood Glucose Monitoring System is intended for use outside the body (in vitro diagnostic use) at home. It is used for quantitative measurement of glucose level in fresh capillary whole blood samples (from the finger, the palm, the forearm, the upper arm, the calf and the thigh). The alternative site testing can be only used during steady-state blood glucose monitoring. The Advocate® Redi-Code<sup>+</sup> blood glucose monitoring system is intended for use by a single person and should not be shared. In addition, it is intended for use at home as an aid in monitoring the effectiveness of diabetes control program. It should not be used for the diagnosis or screening of diabetes, nor for the testing of neonates.

The Advocate® Redi-Code<sup>+</sup> Blood Glucose Monitoring System consists of the Advocate® Redi-Code<sup>+</sup> Blood Glucose meter and the Advocate® Redi-Code<sup>+</sup> Blood Glucose test strips. The Advocate® Redi-Code<sup>+</sup> Blood Glucose meter is used only with Advocate® Redi-Code<sup>+</sup> Blood Glucose test strips to quantitatively measure glucose in fresh capillary whole blood samples drawn from finger tips, the palm, the forearm, the upper arm, the calf and the thigh.

The Advocate® Redi-Code<sup>+</sup> control solutions are for use with the Advocate® Redi-Code<sup>+</sup> Blood Glucose Monitoring system as a quality control check to verify the accuracy of blood glucose test results.

### ***Test Principle***

Advocate® Redi-Code<sup>+</sup> Blood Glucose Test is based on the quantitative measurement of glucose in whole blood. The blood sample is drawn into the tip of the test strip before testing. Glucose in the sample mixes with special chemicals on the test strip and a small electrical current is produced. The current is measured by the meter and the corresponding blood glucose level is displayed. The strength of the current produced by the reaction depends on the amount of glucose in the blood sample.

## **Cautions**

- Only use Advocate® Redi-Code<sup>+</sup> Test Strips with Advocate® Redi-Code<sup>+</sup> Meter. Other test strips may produce inaccurate results.
- Perform all quality control checks recommended in this Owner's Manual.
- As with all small parts in this system, these parts should be kept away from small children who still put things in their mouths. If they are swallowed, promptly see a healthcare professional.
- Read all instructions thoroughly and practice the test before using the product to test your blood glucose. Do all quality control checks as directed and consult with a diabetes healthcare professional.

## Introduction of Your System

### System Contents

#### KIT CONTENTS


- Meter


- Owner's Manual

- Log book

- Test Strips

- Quick Reference Guide


- Control Solution

- Lancet

- Lancing Device with a Clear Cap

#### Notes:

- If this system has been opened prior to use, please return it to the place of purchase.
- For information on how to obtain accessories including test strips, lancing device, lancets and control solutions. Please contact your distributor at 1-866-373-2824 or place of purchase for availability.

## Meter components

### Test Strip Port


Insert test strip. The meter turns on when you insert the test strip.

### Display

Display results and other important information.


### Main Button

Press to enter memory mode to review stored results or adjust setting.


### Battery Cover

Remove cover and replace battery.

### Set Button

Press to enter setting mode and change the date, time, alarm settings and memory deletion.

# Introduction of Your System

## Meter Display


### Year/Date

### Memory

Appears when test results stored in the memory.

### Day Average

Appears when viewing averages (7-, 14-, 28-, 60-day).

### Alarm Marker

Indicate an alarm has been set.

### Time

### Pre-Meal Marker/ Post-Meal Marker

### Control Solution Marker

Indicates the result is performed under control solution mode.

### Battery

Indicates the batteries are low and need to be changed.

### Unit of Measurement

### Test Strip Marker

Appears when the meter is ready for inserting a test strip.

### Blood Drop

Appears when the meter is ready for applying sample.


### KETONE?

Appears when a blood glucose test falls above 240 mg/dL

### Test Result Area

Test results are displayed here .

## Test Strip


## Meter Set-Up


### Install Batteries

Make sure that the meter is turned off before replacing the batteries.


1. Press the buckle on battery cover and lift up to remove cover.

2. Remove the old battery and replace with two 1.5V AAA alkaline batteries.

3. Close the battery cover.


## Notes

- Batteries should be replaced when the power of your meter is getting low by displaying battery symbol or the meter briefly shows “E-2” with battery symbol.

- After replace new batteries, the display check will appear briefly. That means the batteries are installed successfully.

- Replacing the batteries does not affect the meter's memory; however, the meter's setting (date, time and alarms) may need to be re-set.
- Batteries might leak chemicals if not used for a long time. Remove the batteries if you are not going to use this system for an extended period (i.e., 3 months or more).

## Set the Meter

Your meter allows you to set the year, month, date, time, alarm as well as delete memory.

- **Enter Setting Mode**

**Start with the meter off.** Press Set button (the back of the meter).

- **Exit Setting Mode**


Any time you want to exit setting mode, press and hold M button for 2 seconds until the meter turns off.

- **Move Setting Faster**

Hold M button down to scroll forward quickly until the number you desire appears.

### 1. Set the Year

Press Set button to enter Setting mode. The Year will be flashing. Press and release M button until the correct year appears. Press Set button to move to Month setting.


### 2. Set the Month

With the month flashing, press and release M button until the correct month appears. Press Set button to move to Day setting.


### 3. Set the Day

With the day flashing, press and release M button until the correct day appears. Press Set button to move to Hour setting.

### 4. Set the Hour

With the hour flashing, press and release M button until the correct hour appears. Press Set button to move to Minutes setting.

### 5. Set the Minutes


With the minutes flashing, press and release M button until the correct minute appears. Press Set button to move to Alarm setting.

## 6. Set the Alarm

The meter provides six alarms (1-6). You may set up any or all of them.

With the first alarm flashing (the default is OFF), press M button to switch "OFF" to "ON" and go to **step 6-1** to see how to set the alarm's time. Or with "OFF" displaying, keep pressing Set button to move to Memory Deletion.

### 6-1. Enter the Alarm Setting

01 (Alarm 1) and  with flashing "On" displays, press Set button to enter the Hour setting.


### 6-2. Set the Alarm's Hour

With the hour flashing, press and release M button until the correct hour appears. Press Set button to move to Minutes setting.


## 6-3. Set the Alarm's Minutes

With the minutes flashing, press and release M button to select 15 minutes, 30 minutes, or 45 minutes. Press Set button to finalize Alarm 1 setting and then move to Alarm 2 setting. If you want to leave alarm setting, when "OFF" displays with alarm marker, press Set button to move to Memory Deletion.

## 6-4. Set Other Alarms (2-6)

Repeat steps 6-1, 6-2 and 6-3 to set Alarm 2, Alarm 3, Alarm 4, Alarm 5 and Alarm 6. After finalizing Alarm 6, press Set button to move to Memory Deletion.

## 7. Delete Memory


With "M" symbol and flashing "dEL" appears, if you **DO NOT** want to delete memory, press Set button to skip this step. If you want to delete ALL memory, press M button for 2 seconds to clear ALL memory.

## Meter Set-Up

With “CLr” displaying, that means ALL memory is deleted.

After pressing Set button or deleting memory, set-up is all completed. The meter will show “OFF” and then turn off.


English

### Notes

- The date, time and alarm can **ONLY be changed** in the setting mode. Therefore, when performing a glucose testing, those parameters are not possible to be changed.
- The meter displays 7-, 14-, 28- and 60-day averages which you can access from the meter memory. These averages are calculated from results obtained during the 7-, 14-, 28- and 60-days preceding the current date and time settings. When the date and time are changed, the 7-, 14-, 28- and 60-day averages will be re-calculated and updated automatically based on the new date.
- The date and time might need to be adjusted due to daylight saving
- The setting mode needs to be adjusted due to daylight saving and after you replace the battery.
- If no button is pressed for 2 minutes, the meter will turn off automatically.


## Use the Alarm

When the alarm sounds, the meter will automatically turn on. The alarm will keep sounding for 30 seconds. You can perform a test or press M button to turn off the meter.

## Perform a Test

Insert a test strip into the test strip port and follow the Blood Glucose Test instructions (page 25).

## Turn off the Alarm

Press M button to turn off the meter or leave the alarm sounding for 30 seconds and then the meter will turn off.

## Cancel the Alarm

If you want to cancel the alarm function, press Set button when the meter is off. Follow "Set the Meter" instruction to cancel the alarm function (page 16).

## Quality Control

Advocate® Redi-Code<sup>+</sup> Level 1, Level 2 and Level 3 Control Solutions are able to check the meter and the test strips are working together as a system. It is very important that you do this simple check routinely to make sure you results are accurate.

### *When to check your system with control solutions*

- Before using the meter for this first time.
- When you begin using a new vial of test strips.
- Whenever you suspect that the meter or test strips are not working properly.
- After cleaning and disinfection of the meter.
- At least once a week.
- To check your technique.
- When your blood glucose test results are not consistent with how you feel, or when you think your results are wrong.
- When your test strips are exposed to extreme environmental conditions (below 39.2°F, or above 104°F, or humidity levels above 85%).
- If you drop the meter.

### *Important Control Solution Test Information*

- **The control solution range printed on the test strip vial is for Advocate® Redi-Code<sup>+</sup> Control Solution only.** It is used to evaluate meter and test strip performance. It is not recommended range for your blood glucose level.
- Use only **Advocate® Redi-Code<sup>+</sup>** Control Solutions. Check the expiration date on the control solution vial. Do not use if expired.
- Control solution, meter, and test strips should come to room temperature between 68°F and 77°F before testing.
- Shake the vial before use. Discard the first drop of control solution. After squeezing, wipe off the dispenser tip to avoid contamination. The above tips ensure you to get a good sample and an accurate result.
- Due to the potential biohazardous concern, you should wash hands


thoroughly with soap and water after handling the meter, lancing device, or test strips.

- Use for only 90 days from the date that you first open it. Record the discard date (date opened plus 90 days) on the control solution vial. Discard after 90 days.
- Store the control solution tightly closed at temperatures between 35°F-86°F. Do not freeze.

## Run a Control Solution Test

### 1. Insert a Test Strip


Fully insert a test strip to the meter with the contact bars end first. Push it in until it will go no further. The meter will turn on and the display check will appear briefly.

### 2. Press M Button for CTL mode


When the flashing apply blood symbol appears on the display, press M button thrice until a control bottle with "CTL" appears.

**Control results marked with a control bottle and "CTL" will NOT be stored in the meter's memory.**


### 3. Obtain Control Solution

Shake the control solution vial well. Remove the cap and squeeze the vial to discard the first drop. Then wipe off the dispenser tip to prevent contamination. Squeeze the vial again to produce another drop and **place this drop on a clean surface.**


## 4. Apply Control Solution

Hold the meter to touch the drop until the meter beeps and begins to count down.


To avoid contaminating the control solution with the content of the test strip, **DO NOT** directly apply control solution into a strip.

## 5. Read the Control Solution Result

After 5 seconds, the control solution test result appears. Compare the result with the range printed on the test strip vial. The result should fall within this range.

Remove the used strip, the meter will automatically turn off

## Cautions

- The control solution range printed on the test strip vial is for ADVOCATE® Redi-Code+ Control Solution only. It is not a recommended range for your blood glucose test results.
- Every time you perform a control solution test, you must enter into the "CTL" mode so that the result will not be stored in the meter memory. Failure to do so will confuse the blood glucose result with the control solution result in meter's memory.

***If the result is not in the acceptable range, it may be caused by***

- Performing the test improperly.
- Not shaking the control solution vial well enough.
- Using expired or contaminated control solution and test strips.
- Using the first drop of control solution.
- Using a damaged meter.
- Control solution is stored at environment out of the specified temperature range, which is from 35°F to 86°F.

If you continue to get control solution test results that fall outside the range, the system may be not working properly. **DO NOT** use the system to test your blood. **Contact your distributor at 1-866-373-2824 for help.**

## Blood Glucose Test

Make sure you have all items needed to test before you perform your blood glucose.

- Blood Glucose Meter
- Test Strips
- Lancing Device
- Sterile Lancets

### Warnings


To reduce the chance of infection

- Make sure to wash and dry the puncture site before starting.
- Never share a lancet or the lancing device with anyone to avoid the risk of blood borne pathogen transmission.
- Always use a new, sterile lancet. Lancets are for single use only.
- Keep the meter, the lancets and the lancing device clean before using.
- The user should wash hands thoroughly with soap and water after handling the meter, lancing device, or test strips.
- If the meter is being operated by a second person who is providing testing assistance to the user, the meter and lancing device should be decontaminated prior to use by the second person. Please refer to the section "Cleaning and Disinfection" for details in this manual.

### Prepare for the Test


#### 1. Insert a Lancet

Remove the cap and then insert a lancet into the holder and push in firmly.


#### 2. Remove the Protective Disk

Twist the lancet disk cover until it separates from the lancet and leave the cover aside for later use.


### 3. Replace the Lancing Device Cap

Replace the lancing device cap by twisting it back on until it is snug.


### Use the Clear Cap for Alternate Site Testing (AST)

Twist on the clear cap until it is snug.


### 4. Set the Lancet Depth

Twist the adjustable tip in either direction until the number lines up with the arrow. The adjustable tip offers 6 levels of skin penetration. Start at 3 or 4. For tougher skin, twist to a higher number.

### 5. Cock the Lancing Device

Slide the ejection/cocking control back until it clicks. If it does not click, it may have been cocked when the lancet was inserted.

**The lancing device is now ready for use.**


# Blood Glucose Test

## Run a Blood Glucose Test

### 1. Insert a Test Strip

Fully insert a test strip to the meter with the contact bars end first. Push it in until it will go no further.

The meter will turn on and the display check will appear briefly.


### 2. Choose Pre-Meal/Post-Meal/General Mode

You can choose three modes by pressing M button, if you miss the desired mode, you can continue pressing M button to go back.

#### Pre-Meal Mode

When the flashing apply blood symbol appears on the display, pre-meal marker (☞) displays firstly. Go to Step 3 for obtaining a drop of blood


#### Post-Meal Mode

When the flashing apply blood symbol appears on the display, press M button once and post-meal marker (☾) displays. Go to Step 3 for obtaining a drop of blood

#### General Mode

When the flashing apply blood symbol appears on the display, press M button twice and general mode appears (no marker displays). Go to Step 3 for obtaining a drop of blood.


### 3. Obtain a Drop of Blood

Before puncturing, clean the puncture site with 70% alcohol cotton and let it air-dry.


#### Fingertip

Hold the lancing device firmly against the side of your fingertip and press the release button. After hearing a click, it indicates that the puncture is complete.

#### Sampling for an alternate site (AST) other than fingertip

Press and hold the lancing device against your palm, forearm, upper arm, calf or thigh for a few seconds, then press the release button. Keep holding the lancing device and cap against your skin until a round drop of blood appears.


#### Warnings

- The risk of transmission of disease such as viral hepatitis from shared use of lancing device has been increasingly reported by the FDA and CDC<sup>\*1,2</sup>. To prevent the transmission of disease and infection, the lancing device is intended only for a single user and should not be shared.
- Squeeze the punctured area gently to obtain blood and be careful NOT to smear the blood sample.
- Choose a different spot each time you test. Repeated punctures in the same spot may cause soreness and calluses.
- Before you decide to sample for alternate sites, please refer to "Alternate Site Testing" for limitations or consult your healthcare professional.
- It is recommended that you discard the first drop of blood as it might contain tissue fluid and serum, which may affect the test result.

# Blood Glucose Test

## 4. Apply Blood to the Test Strip

Apply the tiny blood sample to the test strip. The minimum blood sample size is 1.1µl (●). Blood will be drawn into the strip. Keep holding the drop of blood to the hole of the test strip until the confirmation window is full and the meter begins to count down. **Please note that if the confirmation window is not filled with blood, the test result may inaccurate.**


## 5. Get Your Result

After 5 seconds, your blood glucose test result appears. This result will be automatically stored in the meter's memory.

## 6. Dispose the Used Lancet Remove the lancet carefully to prevent injury.

Take the lancet out from the lancing device. Place the disk on a hard surface and push the exposed tip into the disk cover.

## Warnings

- You should wash hands thoroughly with soap and water after handling the meter, lancing device, or test strips.
- Please refer to the section of Cleaning and Disinfection for detailed information.
- If no action is done for 2 minutes, the meter will automatically turn off. You must remove and reinsert the test strip to restart the test procedure.
- The blood should be completely filled the confirmation window before the meter begins to count down. If you find that the confirmation window is not filled with blood when the meter is counting, **NEVER** try to add more blood to the test strip. **Discard the test strip and retest with a new one.**
- Any trouble filling the confirmation window, **contact your distributor at 1-866-373-2824 for help.**
- The used lancet and the used test strip may be potentially biohazard. Please discard it carefully according to your healthcare professional's instructions.

## Expected Blood Glucose Values

Many factors will affect blood glucose values such as food intake, medication dosages, health, stress, or exercise. The expected blood glucose values appropriate for you should be consulted and determined by your healthcare professional.

People without diabetic blood glucose levels are expected to be in the following ranges<sup>\*6</sup>:

<i>Time</i>	<i>Range (mg/dL)</i>
Fasting and before meals	70-130
<b>1-2 hours after meals</b>	Less than 180

## References

- 6.American Diabetes Association (2010). Standards of Medical Care in Diabetes. Diabetes Care, 33 (Supplement 1): S11–S61.

### Limitations

- a)** The test strips are used for fresh capillary whole blood samples.
- b)** DO NOT use neonate blood sample.
- c)** Not to be used for diagnosis or screening of diabetes.
- d)** Alternative site testing with this system can be used only during steady-state blood glucose conditions.
- e)** Measurements from alternative site testing should never be used to calibrate a continuous glucose monitor (CGM) or entered into insulin dose calculators for insulin dosing recommendations.
- f)** Extreme humidity may affect the results. A relative humidity greater than 85% may cause incorrect results.
- g)** The system should be used at a temperature between 50°F - 104°F (10°C - 40°C). Outside this range, the system halts.
- h)** DO NOT reuse the test strips. The test strips are for single use only.
- i)** Red blood cell count (Hematocrit); Hematocrit levels between 20% and 60% will not significantly affect the results, Hematocrit levels below 20% may cause higher results, Hematocrit levels above 60% may cause lower results. If you do not know your hematocrit level, please consult with your healthcare professional.
- j)** Altitude up to 10,745 feet above sea level has no effect on readings.

## **Comparison of Your Meter and Laboratory Results**

Test results with this meter are plasma-calibrated. Since the laboratory tests provide plasma-calibrated results, your meter may differ from your laboratory result due to normal variation. Meter results can be affected by factors and conditions that do not affect laboratory results in the same way. To obtain an accurate comparison between meter and laboratory results, follow a few basic guidelines below.

### **Before going to the lab**

- Perform a control solution test to make sure that the meter is working properly.
- Fast for at least eight hours before you perform your blood test.
- Take your meter with you to the lab.

### **While at the lab**

- Only use fresh capillary blood obtained from the fingertip.
- Collect your blood into heparin-containing test tubes and do NOT use other anticoagulants or preservatives.
- Conduct your meter test within 15 minutes of the lab test.

You may still have a variation from the result because blood glucose levels can change in some specific situations that could cause a difference of more than  $\pm 20\%$ .

- You have eaten recently. The blood glucose level from a finger stick can be up to 70 mg/dL higher than blood drawn from a vein (venous sample) used for a lab test<sup>\*7</sup>.
- If your hematocrit (the amount of red blood cells in the blood) is high or low.
- You are suffering from severe dehydration (the loss of body fluid).

You may see the package insert of test strip for accuracy and precision data and for important information on limitations.

## ***References***

7. Sacks, D.B.: "Carbohydrates." Burtis, C.A., and Ashwood, E.R.( ed.), Tietz Textbook of Clinical Chemistry. Philadelphia: W.B. Saunders Company (1994), 959.

## Meter Memory

The meter stores the last 400 blood glucose test results along with the date and time. It also provides you with 7-, 14-, 28- and 60-day averages of your pre-meal, post-meal and general blood glucose test results, respectively. You can begin to review the results stored in meter's memory by following the next steps.

### Enter Memory Mode

**Start with the meter off.** Press M button until a test strip marker appears.


### Exit Memory Mode

Any time you want to exit memory mode, press and hold M button for 2 seconds until the meter turns off.

### 7-,14-,28-,60-Day Average Results

#### 1. Enter Memory Mode

When the meter is off, press and release M button. The meter will show flashing strip symbol, it means you can test your blood glucose right now. If you do not want to test a blood, press M button to review 7-day average result for pre-meal mode.


#### 2. 7-Day Average Pre-Meal Result

The 7-day average pre-meal result then appears after press and release M button.

**07 DAY AVG** = your 7-day average

**30** = indicates the number of results included in your average

 = pre-meal marker

## 3. 7-Day Average Post-Meal Result

Press and release M button, then 7-day average post-meal result appears.

**07 DAY AVG** = your 7-day average

**30** = indicates the number of results included in your average

☾ = post-meal marker

## 4. 7-Day Average General Result

Press and release M button, then 7-day average general result appears.


**07 DAY AVG** = your 7-day average

**30** = indicates the number of results included in your average


## 5. 14-Day Average Results for Pre-Meal, Post-Meal and General Mode

Press M button again, 14-day average pre-meal result appears. Like Step 2 to 3, 14-day average post-meal result and general result will appear in order.


## 6. 28-Day Average Results for Pre-Meal, Post-Meal and General Mode

After 14-day average general result, press M button, 28-day average pre-meal result appears. Like Step 2 to 3, 28-day average post-meal result and general result will appear in order.

## 7. 60-Day Average Results for Pre-Meal, Post-Meal and General Mode

After 28-day average general result, press M button, 60-day average pre-meal result appears. Like Step 3 to 4, 60-day average post-meal result and general result will appear in order.

Continue to press the M button to view stored test results. The most recent result will display first.

## Test Results

### 1. View Test Results

After viewing 60-day average result, press M button again. The most recent result with date and time will be shown on the display.

### 2. View Second Most Recent Result

Press the M button once by once to recall the test results stored in the meter consecutively.



## 3.View the Last Result

After viewing the last result, press M button again and the meter will be turned off.

### Note

- The control solution results are NOT stored in the memory. The list of past results and the result average are for blood glucose results only.
- The stored results marked with a pre-meal marker (☞)/post-meal marker (☛) depend on what you chose in your test procedure. If you chose general mode in your test procedure, the stored results without markers will display.
- If no button is pressed for 2 minutes, the meter will show “OFF” and turn off automatically.
- If no results are stored in the meter's memory, “---” will be shown on the LCD.


## Your System Maintenance

### Cleaning and Disinfection

The risk of transmission of disease such as viral hepatitis from shared use of the meter and lancing device has been increasingly reported by the FDA and CDC<sup>\*1,2</sup>. The cleaning procedure is to thoroughly remove all debris, blood and other body fluids from surfaces and the disinfection procedure is to destroy or prevent growth of organisms capable of infection. To prevent the transmission of disease and infection, **the Advocate<sup>®</sup> Redi-Code<sup>+</sup> Blood Glucose Meter and Lancing Device should be cleaned and disinfected when the following situations happen.**

- **Whenever the meter and lancing device are visibly dirty.**
  - **At least once per week.**
  - **Before the meter and lancing device are being operated by a second person who is providing testing assistance to the user.** Make sure you have all items before you clean and disinfect the meter and lancing device.
  - A new pair of clean gloves, which are available in pharmacies near your location.
  - A commercially available disinfectant towelette, CaviWipesXL (manufactured by Metrix with EPA registration number 46781-8). Please contact Metrix Research Corp. at 1-800-841-1428 for availability.
- \* The CaviWipesXL has been shown to be safe for use with Advocate<sup>®</sup> Redi-Code<sup>+</sup> Blood Glucose meter and lancing device, but any other disinfectant product with the EPA registration number 46781-8 may be used on this device.


## Area to be cleaned and disinfected for the meter

## Possibility of contact with blood for all parts of meters

Part	Name of the part	Possibility of contact with blood
A	Display	Low
B	Strip insert port	High
C	Meter housing	High
D	Battery cover	Low
E	Main button	Low

## Area to be cleaned and disinfected for the lancing device

## Possibility of contact with blood for all parts of the lancing device

Part	Name of the part	Possibility of contact with blood
A	Lancing device cap	High
B	Clear cap	High
C	Main body of the device	Low
D	Adjustable tip	Low
E	Lancet ejector	Low


Do	Do not
<ul style="list-style-type: none"><li>* Make sure the meter is turned off during cleaning and disinfection.</li><li>* Keep the test strip vial tightly closed when performing the cleaning and disinfection procedures because the fumes from the disinfectant may affect the performance of the strip.</li><li>* After cleaning or disinfection, please perform the physical appearance and performance check of devices.</li></ul>	<ul style="list-style-type: none"><li>* Get any moisture in the test strip slot.</li><li>* Spray any cleaning solution directly onto the meter.</li><li>* Put the meter under water or liquid.</li><li>* Pour liquid into the meter.</li></ul>

### 1. Cleaning Process

**1-1.** Put on disposable protective gloves and take disinfectant CaviWipesXL from the container.

**1-2.** Wipe down the body of the meter and lancing device. After air-drying, go to Step 2, Disinfection for the meter and lancing device.

**1-3.** Discard the used disinfectant wipe in trash.


## 2. Disinfection Process

\* A pre-cleaning step (using the CaviWipesXL) is required prior to each disinfection step.

**2-1.** Take a new disinfectant CaviWipesXL from the container.

**2-2.** Wipe down the body of the meter and lancing device enough to thoroughly wet the surface. Repeat this step for six times (3 passes horizontal and 3 passes vertical as shown in the figure) to complete the disinfection procedure. Allow the devices to remain wet for a 2-minute contact time. Similarly clean and disinfect the lancing device, including the clear cap used for alternate site testing. Do disinfection at least once per week .

**2-3. Discard the used disinfectant wipe in trash.**

**2-4. The user should wash hands thoroughly with soap and water after handling the meter, lancing device, or test strips.**

## Your System Maintenance

After each cleaning or disinfecting process, you should perform the physical appearance check and performance check of the meter and lancing device for items shown in the tables below.

### Physical appearance check-after each cleaning and disinfection for the meter

Check item	Acceptance Standard
Losing display transparency	Accept the result if none of the items listed on the left appears after each cleaning or disinfection process
Losing legibility of any marking	
Polymer crazing	
Cracking	
Swelling	
Dissolving	
Softening	
Becoming brittle	
Deformation	
Cloudy	

### Physical appearance check-after each cleaning and disinfection for the meter

Steps	Check items	Accept result
1. Insert test strip	Is the meter powered on?	Accept result if the meter is powered on

2.Press M Button for CTL mode	Is “CTL” showing on the display?	Accept result if the “CTL” is showing on the display
3.Apply control solution to obtain the control solution test result	Is the control test result within the control solution range specified on the strip vial?	Accept result if the control test result is within the control solution range specified on the strip vial

**Physical appearance check- after each cleaning or disinfection for the lancing device, including normal and AST cap**

Check item	Acceptance Standard
Polymer crazing	Accept the result if none of the items listed on the left column appears after each cleaning or disinfection process
Cracking	
Swelling	
Dissolving	
Softening	
Becoming brittle	
Deformation	
Cloudy	


## Performance check-after each cleaning or disinfection for the lancing device, including normal and AST cap

Steps	Check items	Accept result
1.Capping of normal and AST cap	Could the normal and AST cap fit perfectly with the main body of the device?	Accept result if both caps fit perfectly with the main body of the device
2.Cock the device and press on the trigger	Does the device work after pressing the trigger?	Accept result if the device works
3.Install a lancet to the lancet holder and push in firmly	Could the lancet holder hold the lancet firmly?	Accept result if the holder could hold the lancet firmly

In case the result is not acceptable, please contact customer care service, with the contact information listed on the last page of the manual, for replacement with a new device.

## Note

- Special attention should be given during the cleaning and disinfecting processes- not to squeeze the disinfectant, which might cause drop of the virucide to get into the test strip insertion hole of the meter and thus cause the meter to malfunction.

Virucide cannot drop into the strip insert hole.

- You should wash hands thoroughly with soap and water after handling the meter, lancing device, or test strips.
- Any trouble performing cleaning and disinfection procedure on the meter and lancing device, **please contact your distributor at 1-866-373-2824 for technical assistance.**
- Based on the robustness study performed by the manufacturer, the meter has passed all the appearance and performance checks after a total of 15,000 cleaning and disinfection cycles, and the lancing device has passed all the appearance and performance checks after a total of 5,000 cleaning and disinfection cycles. The meter must be replaced with a new meter at the end of 15,000 cleaning and disinfection cycles or 5 years, whichever comes first. The lancing device must be replaced with a new lancing device at the end of 5,000 cleaning and disinfection cycles or 1 year, whichever comes first.

## Storage

- Always store or transport your meter, test strips and control solution in the carrying case after each use.
- Store your meter, test strips and control solution in a cool and dry place. Keep away from direct sunlight and heat.
- Avoid dropping and strong impact on your system.
- After removing a test strip from the vial or using control solution, immediately replace the vial cap and close it tightly to avoid contamination or damage.
- Do not use test strips and control solution after expiration date. Write the discard date on their vials after first opening.
- Keep the strip vial away from children since the cap and the test strip may be a choking hazard. If swallowed, promptly see a healthcare professional for help.





## Messages and Troubleshooting

This system displays messages when there are problems with the test strip, with the meter, or when your blood glucose levels are lower than 20 mg/dL or higher than 600 mg/dL. These messages help to identify certain problems but do not appear in all cases when a problem has occurred. Improper use may cause an inaccurate result without producing an error message. In the event of a problem, refer to the information under "What to do".

Message	What it means	What to do
	Your blood glucose result is lower than 20 mg/dL (below measurement limit).	You should seek immediate medical assistance.
	Your blood glucose result is higher than 600 mg/dL (above measurement limit).	You should seek immediate medical assistance.
	Your blood glucose result is equal to or higher than 240 mg/dL, <b>KETONE?</b> appears with your result.	Please seek immediate medical assistance.

# Messages and Troubleshooting

Message	What it means	What to do
 <p>The display shows "E-1" and a battery icon with a horizontal bar inside, indicating a low battery level.</p>	<p>The test strip cannot be identified. A previously used test strip may have been inserted.</p>	<p>Use a new test strip to perform test strip. Be sure you use ADVOCATE<sup>®</sup> Redi-Code<sup>+</sup> test strips.</p>
 <p>The display shows "E-1" and a blood sample icon with a horizontal bar inside, indicating an insufficient blood sample.</p>	<p>Insufficient blood sample or the test strip is not completely filled with blood sample.</p>	<p>Remove the test strip and retest using a new test strip. Be sure you apply sufficient blood sample to fill test strip.</p>
 <p>The display shows "E-2" and a battery icon with a horizontal bar inside, indicating dead batteries.</p>	<p>Batteries are dead.</p>	<p>Replace the batteries immediately.</p>
 <p>The display shows a battery icon with a horizontal bar inside, the date "03-10", the time "06:38 AM", the test result "137", and the unit "mg/dL".</p>	<p>When battery symbol displays with the test result, battery power level is low.</p>	<p>Replace the batteries as soon as possible.</p>

Message	What it means	What to do
	Operating temperature is out of range.	Move to an area between 50°F and 104°F, let the meter adjust to this temperature for 20 minutes before performing a test.
	No results are stored in the meter's memory.	No action is required.
	No average results are stored in the meter's memory.	No action is required.
	The alarm sounds to remind you of performing a test.	Insert a new test to perform a test or press M button to turn off the meter.
<b>Meter does not display a message after inserting a test strip.</b>	✓ Possible dead battery	Replace the battery.
	✓ Battery incorrectly installed or missing.	Check that the battery is correctly installed.
	✓ Test strip inserted upside down or incompletely.	Insert the test strip correctly with the contact bars end.
	✓ Defective meter.	Please call your distributor for help.

# Messages and Troubleshooting

Message	What it means	What to do
<b>Test does not start after applying the sample.</b>	✓ Insufficient blood sample.	Repeat the test using a new test strip with larger volume of blood sample.
	✓ Defective test strip.	Repeat the test with a new test strip.
	✓ Sample applied after automatically shutoff (2 minutes after last user action).	Repeat the test with a new test strip. Apply sample only when blood drop symbol appears on the display.
	✓ Defective meter.	Please call your distributor for help.
<b>Control solution test result is out of range.</b>	✓ Error in performing the test.	Read the instruction thoroughly and repeat the test again.
	✓ Do not shake the control solution vial very well.	Shake the control solution vigorously and repeat the test again.
	✓ Expired or contaminated control solution.	Check the expiration date or the discarded date of the control solution.
	✓ Control solution that is too warm or too cold.	Control solution, meter, and test strips should come to room temperature (68-77°F) before testing.
	✓ Test strip deterioration.	Repeat the test with a new test strip.
	✓ Defective meter.	Please call your distributor for help.
	✓ Control solution is stored at environment out of the specified temperature range, which is from 35°F to 86°F.	Use a new bottle of control solution that has been stored at specified temperature range.

Never try to disassemble the meter in any circumstances. If you encounter any error messages not listed above or if you have followed the actions recommended but the problem remains unsolved, please call your distributor at 1-866-373-2824 for support.

## Specifications

<b>Model No.</b>	BMB-EA001
<b>Dimension</b>	64mmx95mmx29mm
<b>Weight</b>	53g without batteries
<b>Power Source</b>	Two 1.5V AAA alkaline batteries
<b>Battery Life</b>	Over 1,000 tests
<b>Auto Turn-off</b>	2 minutes after idle
<b>Test Sample</b>	Capillary whole blood
<b>Sample Size</b>	1.1 microliter
<b>Test Time</b>	5 seconds
<b>Measurement Range</b>	20 to 600 mg/dL
<b>Meter/Strip Storage Condition</b>	35°F to 86°F, < 85% R.H. (noncondensing)
<b>Operating Condition</b>	50°F to 104°F, 15% / 85% R.H. (noncondensing)
<b>Memory</b>	400 blood glucose tests with date and time
<b>Control Solution Storage Condition</b>	35°F to 86°F

The device has been certified to meet the electrical and safety requirements of: IEC 61010-1, EN 61010-1, IEC 61010-2-101, EN 61010-2-101, IEC 61326-1, EN 61326-1, IEC 61326-2-6, EN 61326-2-6.

## **Statements**

- The equipment complies with IEC 61326-1:2005. Emissions are low and not likely to cause interference in nearby electronic equipment.
- The equipment is tested for immunity to electrostatic discharge as specified in IEC 61000-4-2, at a level of 8kV for Air Discharge and 6kV for Contact Discharge.
- The equipment is tested for immunity to radio frequency interference at the frequency range of 80 MHz to 2.7GHz.

## Warranty

Suncoast, Inc. warrants to the original purchaser that this instrument will be free from defects in materials and workmanship for a period on ONE (1) year from the date of purchase (except as noted below). The warranty is not transferable.

### **This warranty is subject to the following exceptions and limitations**

- 1.This warranty is valid only if the Warranty Registration Card is completed with date of purchase and serial number and sent to Suncoast's office within 14 days commencing from the date of purchase.
- 2.This warranty is limited to replacement due to defects in parts and workmanship.
- 3.This warranty does not apply to the performance of this product that has been damaged due to abuse, accidents, alteration, misuse, neglect, maintenance by someone other than Suncoast's, or failure to operate the instrument in accordance with instructions. Further, Suncoast's assumes no liability for malfunction or damage caused by the use of reagents other than reagents manufactured or recommended by Suncoast.
- 4.Suncoast reserves the right to make changes in design of this instrument without obligation to incorporate such changes into previously manufactured instruments.

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## Referencias

- 1.FDA Public Health Notification: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens: Initial Communication (2010) <http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/ucm224025.htm>
- 2.CDC Clinical Reminder: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens (2010) <http://www.cdc.gov/injectionsafety/Fingerstick-DevicesBGM.html>
- 3.NCCLS, Vol. 25, No. 27, 2005, EP7-A2, Interference testing in clinical chemistry; approved guideline – second edition.
- 4.Tietz Textbook of Clinical Chemistry and Molecular Diagnosis (4th Ed.), Burtis, Ashwood & Bruns (Eds), Elsevier Saunders, 2005; 2301.
- 5.Bioavailability of tolazamide from tablets: Comparison of in vitro and in vivo results, P. G. Welling, R. B. Patel, U. R. Patel, W. R. Gillespie, W. A. Craig, K. S. Albert, 2006.

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